endoBarrier®

The EndoBarrier® Gastrointestinal Liner with Delivery System

Instructions for use

Sterile device for one time use only. Do not reuse or re-sterilize. Read instructions before using.

The EndoBarrier Gastrointestinal Liner System can be purchased only upon order of a physician and is intended for use by qualified medical personnel who have received training by GI Dynamics or its designate.

NOT AVAILABLE FOR SALE OR DISTRIBUTION WITHIN THE UNITED STATES. FOR EXPORT ONLY.

Note: The EndoBarrier® System requires use of a commercial gastroscope and exchange length, 0.035 - 0.038-inch nitinol guidewire.

Product Description

The GI Dynamics' EndoBarrier® Gastrointestinal Liner System is used for the treatment of obese type 2 diabetes with BMI ≥ 30 kg/m², or obese patients with BMI ≥ 30 kg/m² with ≥ 1 comorbidities, or obese patients with BMI > 35 kg/m². It is provided sterile and consists of a gastrointestinal (GI) liner (anchor and liner), preloaded in a catheter that delivers the liner to the proximal intestine. The product is sterilized using ethylene oxide. The GI liner is removed using the EndoBarrier Gastrointestinal Liner Retrieval System, consisting of a Grasper and Retrieval Hoods, which are compatible with standard gastroscopes. The Retrieval System may be supplied either sterile or non-sterile.


These instructions pertain only to the EndoBarrier Gastrointestinal Liner with Delivery System. For information regarding the requirements for EndoBarrier removal, refer to the EndoBarrier® Gastrointestinal Liner Retrieval System *Instructions for Use*.

• Instructions for Use	1
• Gebruiksaanwijzing NL	9
• Mode d'emploi FR	17
• Istruzioni per l'uso IT	25
• Gebrauchsanweisung DE	33
• Instrucciones de uso SP (EU)	41
• Instrucciones de uso SP (LA)	49
• Brugsanvisning DA	57
• Návod k použití CZ	65

EndoBarrier GI Liner

Illustration of the GI liner in the duodenum


The EndoBarrier GI Liner has:


- A nitinol anchor
- An impermeable, fluoropolymer liner with a proximal radiopaque marker
- Two retrieval drawstrings

The anchor attaches the GI liner to the wall of the bulbous duodenum, proximal to the ampulla of Vater. The anchor provides a seal to ensure that chyme passes inside the liner, and barbs on the circumference of the anchor engage the muscularis for fixation. The EndoBarrier liner extends 61 cm into the small bowel from the anchor site.

Delivery System

1. GI Liner Deployment Handle
2. GI Liner/Ball Locking Wire
3. Stiffening Wire Handle
4. Anchor Locking Wire
5. Anchor Deployment Handle
6. Outer Catheter
7. Inner Catheter
8. Capsule

The EndoBarrier GI liner is pre-loaded into and delivered on a custom Delivery Catheter. The Delivery Catheter is approximately three (3) meters long for trans-oral delivery of the EndoBarrier into the jejunum. The catheter is sufficiently flexible to track through the intestine and has a ball on its distal end to facilitate atraumatic tracking through the bowel. The collapsed GI liner is packed into a capsule at the distal end of the outer catheter. Once the EndoBarrier is fully extended in the intestine as seen on fluoroscopy, the distal ball is released and passes through the bowel naturally. The anchor is then deployed in the bulbous duodenum and the catheter is removed.

Retrieval System

The EndoBarrier GI liner is designed to be removed using the EndoBarrier® Gastrointestinal Liner Retrieval System. A separate *Instructions for Use* is available for the retrieval system.

Indications for Use

The GI Dynamics' EndoBarrier® Gastrointestinal Liner System is used for the treatment of obese type 2 diabetes with BMI ≥ 30 kg/m², or obese patients with BMI ≥ 30 kg/m² with ≥ 1 co-morbidities, or obese patients with BMI > 35 kg/m². The GI liner is indicated for a maximum implant duration of 12 months.

Contraindications

The EndoBarrier GI liner is contraindicated in patients with the following:

- Women who are pregnant
- Requirement for prescription anticoagulation therapy
- History of inflammatory bowel disease or condition of the gastrointestinal tract, such as PUD, ulcers, or Crohn's disease
- Pancreatitis
- Uncontrolled GERD
- Known infection at the time of implant
- Symptomatic coronary artery disease or pulmonary dysfunction
- History of coagulopathy
- Bleeding diathesis, upper gastro-intestinal bleeding conditions such as esophageal or gastric varices, or congenital or acquired intestinal telangiectasia
- Congenital or acquired anomalies of the GI tract such as atresias or stenoses
- Previous GI surgery that could affect the ability to place the GI liner or its function
- Iron deficiency and/or iron deficiency anemia
- Active symptoms of kidney stones or known presence of gallstones
- Inability to discontinue use of non-steroidal anti-inflammatory drugs (NSAIDs) during the implant period
- Family or patient history of known diagnosis or pre-existing condition of systemic lupus erythematosus, scleroderma, or other autoimmune connective tissue disorder
- H. pylori positive status (Patients may receive the GI liner if they had a prior history and were successfully treated.)

Warnings

- Only physicians trained in therapeutic endoscopic techniques should use this product.
- Physicians attempting to recover a migrated GI liner should be limited to those who have completed appropriate EndoBarrier retrieval technique training. Devices that have migrated more than 3 cm or have rotated should be removed.
- If the patient exhibits any signs of bleeding that could be related to the device, it is recommended that the GI liner be removed.

- Physicians must be prepared to deal with urgent requirements to remove a GI liner that results in duct or GI liner obstruction.
- Do not withdraw the GI liner through the alimentary tract unless it has been adequately covered by the EndoBarrier Gastrointestinal Retrieval System Retrieval Hood component placed on the distal end of the gastroscope. A liner's anchor that is not adequately covered could damage tissue.
- Remove the gastroscope before removing the outer catheter of the Delivery System; both will not fit in the esophagus. Never allow the outer catheter and capsule to enter the esophagus if the gastroscope is in place as this may lead to perforation.
- Do not attempt to repair or reuse a damaged product, return to manufacturer. Do not reuse or re-sterilize this product. Reuse of this product could cause liner failure or pose a contamination issue.

EndoBarrier Gastrointestinal Liner Precautions

Handling precautions

- Inspect the Delivery System before use to determine that the tracking ball at the tip of the catheter has not separated from the catheter. If it has, gently push the ball back into the capsule.
- Do not modify the system. Modifications could damage the system, resulting in failure to release the GI liner, premature release of the tracking ball, and improper seating of the GI liner in the proximal duodenum.

Pre-procedural precautions

- Patients must receive a proton pump inhibitor (40 mg omeprazole twice a day or equivalent) for three (3) days prior to the procedure and should continue the medication until two (2) weeks after the GI liner is removed.
- Patients should not take anticoagulants (aspirin, heparin, NSAIDs, etc.) for 10 days prior to GI liner placement and for the duration of treatment.
- To reduce the potential for infection, a single 2-gram dose of ceftriaxone (or equivalent) should be administered intravenously 1-2 hours prior to GI liner placement. Individuals with known allergies or hypersensitivity to ceftriaxone, cephalosporins or penicillins should seek an equivalent, long acting, broad spectrum antibiotic. Please see antibiotic package insert for full prescribing information and details.

Procedural precautions

- The ability of the patient to remain still during the procedure must be considered when determining the method of sedation. The type of sedation selected is left to physician judgment. Conscious sedation has not been fully evaluated in a clinical setting.
- A thorough gastroscopic examination of the stomach, pylorus and duodenum should be performed prior to GI liner placement to ensure that the patient's alimentary canal is free of abnormalities which could interfere with the delivery, function and removal of

the liner. This should include observation of the papilla and other accessory ducts. This examination should be made to ensure that they are not located in the bulbous duodenum.

- Placement and removal of the GI liner requires the use of fluoroscopic guidance. Only *water soluble contrast* should be used. Gastrografin or equivalent is recommended. Do not use barium. Caution should be used to protect the patient's reproductive organs from the effects of radiation.
- An incorrectly positioned GI liner may interfere with the bile duct or pylorus. Physicians should remove the incorrectly positioned liner in these patients. Do not attempt to drag or push the liner into position. Additionally, the GI liner should be removed in any patient with clinical and/or chemical evidence of a biliary tract obstruction.
- Do not leave GI liner in place with a crossed tip. Crossed tips should be uncrossed by lifting them with endoscopic forceps.

Post-procedural precautions

It is unknown if certain foods (such as grains and nuts) could interfere with the proper functioning of the gastrointestinal liner. Therefore, patients receiving the GI liner will be placed on a diet regimen similar to that of a Roux-en-Y gastric bypass patient.

Potential Adverse Events

As with all endoscopic and/or implant procedures, serious injury or death can occur.

Potential complications during the placement procedure may include:

- GI tract laceration
- Oropharyngeal perforation
- Esophageal perforation
- Gastric perforation
- Bowel perforation
- Bleeding - in rare cases, significant bleeding that may result in surgery and/or transfusion
- Aspiration
- Infection with or without fever
- Complications associated with endoscopic procedures (i.e., sore throat or deep vein thrombosis)
- Prolonged procedure time
- Headache
- Reaction to antibiotic use
- Anchor across pylorus (misplacement)
- Dizziness

Potential complications during the treatment period may include:

- Small bowel/GI liner obstruction
- GI liner migration or rotation
- Bezoar
- Erosion
- Bleeding - in rare cases, significant bleeding that may result in surgery and/or transfusion
- Vitamin and mineral deficiency

- Dehydration - in rare cases, leading to compromised renal function
- Constipation
- Bloating
- Diarrhea
- Infection with or without fever
- Hepatic abscess
- Pancreatitis (with or without bile duct blockage)
- Gastroenteritis
- Hypoglycemia
- Hyperglycemia
- Flatulence
- GERD
- Esophagitis
- Pseudopolyps
- Nausea/vomiting
- GI pain/cramping
- Peptic ulcer disease
- Duodenitis
- Gastric or bowel perforation
- Local inflammatory tissue reaction
- Flank/back pain
- Alopecia
- Intestinal intussusception
- Headache
- Biliary disease/cholelithiasis
- Hypotension
- Hypertension
- Serum chemistry changes
- Adhesions
- Abdominal pain
- Concurrent Illness
- Weakness (due to diet)
- Patient intolerance

Potential complications during and after GI liner removal may include

- GI tract laceration
- Oropharyngeal perforation
- Esophageal perforation
- Gastric perforation
- Bowel perforation
- Bleeding - in rare cases, significant bleeding that may result in surgery and/or transfusion
- Aspiration
- Adynamic ileus
- Infection with and without fever
- Hepatic abscess
- Pancreatitis (with or without bile duct blockage)
- Biliary disease/cholelithiasis
- Inability to remove GI liner endoscopically, resulting in surgical removal
- Prolonged procedure time
- Complications associated with an endoscopic procedure (i.e., sore throat or deep vein thrombosis)
- Headache
- Alopecia
- Abdominal pain
- Pain—flank/back

- GERD
- Hypoglycemia
- Hyperglycemia
- Concurrent illness
- Diarrhea
- Constipation
- Serum chemistry changes
- Vitamin and mineral deficiency
- GI liner migration/rotation
- Nausea/vomiting
- Weakness

EndoBarrier Gastrointestinal Liner Delivery

How supplied

The EndoBarrier Gastrointestinal liner with delivery system is packaged sterile. The product is sterilized with ethylene oxide.

Equipment Requirements

1. Fluoroscopy unit
2. Video GI endoscopy system, gastroscope with a 2.8 mm working channel and OD of 8.6 – 9.8 mm.
3. Gastrografin, Renografin, or equivalent water soluble contrast media for fluoroscopy, 60 cc. Do not use barium.
4. 60 cc syringes, 20 cc syringes
5. Sterile saline, 200 cc
6. .035 inch super stiff nitinol guidewire
7. EndoBarrier Gastrointestinal Liner with Delivery System (Supplied by GI Dynamics)
8. EndoBarrier Gastrointestinal Liner Retrieval System (Supplied by GI Dynamics)
9. An esophageal overtube 50 cm in length that will fit an 8.6-10 mm endoscope

Patient preparation

1. Patients must eat nothing per mouth for eight (8) hours prior to the procedure.
2. Patients must be started on a proton pump inhibitor (40 mg omeprazole twice a day or equivalent) three (3) days before the procedure and should continue the medication until two (2) weeks after the GI liner is removed.
3. Antispasmodics may be used.
4. A 2-gram dose of ceftriaxone (or equivalent) should be administered intravenously 1-2 hours prior to GI liner placement.

Sedation

The physician should determine if general anesthesia or conscious sedation is appropriate for the patient.

Procedural Steps

1. Lay the patient on the left side.
2. Advance the gastroscope into the duodenum and perform a surveillance exam.
3. Advance the super-stiff guidewire into the duodenum through the gastroscope working channel. Advance and loop the tip of the wire in the duodenum.

4. Remove the gastroscope over the wire while maintaining the wire position. Use fluoroscopy to ensure the loop of wire remains fixed in the duodenum. Minimize the length of wire through the stomach.
5. While maintaining guidewire position, advance the capsule of the catheter over the guidewire and into the duodenum. If needed, pull back slightly on the guidewire to lift the capsule into the pylorus. The gastroscope may be used to help place the capsule into the pylorus.
6. Remove the guidewire. If the guidewire is not easily removed, pull the capsule out of the pylorus into the stomach and try again. If significant resistance is still met, remove the gastroscope and then remove the entire catheter system and start over.
7. **Slowly** advance the inner catheter, by pressing the button on the handle (#1) and sliding the handle piece forward. Release the button to slide the handle piece back and repeat. Advance the inner catheter while observing with fluoroscopy until the catheter is fully delivered as indicated by the distal most reference marker on the inner shaft. If resistance is met, apply forward pressure and wait for peristalsis to advance the catheter.

Notes:

- **Advance the inner catheter slowly. The catheter may prolapse if advanced too quickly.**
 - **The GI liner is fully delivered when the proximal radiopaque marker advances approximately 2-3 cm distal to the capsule.**
8. Retract the inner catheter locking wire (#2) 10 cm to release the distal ball and the liner.
 9. Advance the stiffening wire (#3) to push the ball off the distal end of the catheter. Use fluoroscopy to confirm that the ball is released. Once the ball is off, pull the stiffening wire back (#3) to its stop position.
 10. Reintroduce the gastroscope into the stomach. Ensure that the delivery capsule is fully in the bulb.

Precaution: If the capsule does not fit entirely within the bulbous duodenum, remove the gastroscope then the entire system. This is an indication that the GI liner may not fit well and should not be placed.

11. Retract the anchor locking wire (#4) 10 cm to unlock the anchor.
12. Under endoscopic visualization, position the capsule so that the solid black capsule marker is aligned with the proximal side of the pylorus. Advance the anchor plunger (#5) to deploy the anchor partially out of the capsule, with the barbs remaining in the capsule. Verify the position of capsule, then deploy the anchor by pushing the anchor deployment handle (#5).

Note: Once the anchor is deployed, continue removing the capsule from the bulb through the pylorus while holding the inner catheter steady, until the capsule is in the stomach.

Warning: Do not permit the capsule to advance into the esophagus. Maintain the capsule in the stomach distal to the gastroscope at all times until the gastroscope is removed.

13. Remove the stiffening wire (#3) from the inner catheter.
14. Inject roughly 60 cc of saline or 20% Gastrografin solution through the stiffening wire port of the inner catheter to fill the GI liner with fluid. Then inject

- approximately 60 cc of air through the same port to further separate the inner catheter from the liner.
15. Remove the inner catheter from the outer catheter while watching fluoroscopy to determine if the implant is bunching backward.
 16. Remove the gastroscope from the patient.
 17. Remove the outer catheter from the patient.
 18. Re-introduce the gastroscope into the patient to examine the GI liner in the duodenal bulb.
 19. Flush 60 cc of 20% Gastrografin solution through the working channel of the gastroscope into the GI liner to confirm patency of the liner and also to stretch out the sleeve section of the GI liner.

Note: If the GI liner does not appear to be fully expanded, remove it and insert another liner.

20. Remove the gastroscope.

Notes and Precautions:

- **If anchor tips are crossed, use an endoscopic forceps to lift one leg off the other.**
- **If the anchor is incorrectly positioned, remove the GI liner; do not try to reposition. Remove the gastroscope and then remove the GI liner using the EndoBarrier Gastrointestinal Retrieval System following the Retrieval System *Instructions for Use*.**
- **If the retrieval drawstring is wrapped around the anchor or barbs, remove the GI liner using the EndoBarrier Gastrointestinal Retrieval System following the Retrieval System *Instructions for Use*.**

MR Conditional

The EndoBarrier GI liner has been determined to be MR-conditional. Non-clinical testing has demonstrated that the GI liner is MR-conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 720-Gauss/cm or less

MRI-related Heating

In non-clinical testing, the EndoBarrier GI liner produced the following temperature rise during MRI performed for 15 min of scanning (i.e., per pulse sequence) in the 3-Tesla (3-Tesla/128-MHz, Excite, HDx, Software 14X.M5, General Electric Healthcare, Milwaukee, WI).

Highest temperature change +1.5°C

Therefore, the MRI-related heating experiments for the EndoBarrier GI liner at 3-Tesla using a transmit/receive RF body coil at an MR system reported whole body averaged SAR of 2.9 -W/kg (i.e., associated with a calorimetry measured whole body averaged value of 2.7-W/kg) indicated that the greatest amount of heating that occurred in association with these specific conditions was equal to or less than +1.5°C.

Artifact Information

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the EndoBarrier GI liner. Therefore, optimization of MR imaging parameters to compensate for the presence of the liner may be necessary.

Pulse Sequence	T1-SE	T1-SE	GRE	GRE
Signal Void Size	3.555 mm ²	3.374 mm ²	7.399 mm ²	5.582 mm ²
Plane Orientation	Parallel	Perpendicular	Parallel	Perpendicular

EXPLANATION OF LABEL SYMBOLS

	Consult <i>Instructions for Use</i>
	Sterilized using Ethylene Oxide
	Do not reuse
 YYYY-MM	Use by date or sterility expiration date
	Contains one system
	Do not use if package is damaged
	Catalog number
	Manufacturer's lot number
	Manufacturer
	Authorized Representative in the European Community
	MR-conditional

Warranty

The handling, storage and use of this device relating to the patient, diagnosis treatment, surgical procedures, and other matters beyond GI Dynamics' control directly affect the device and the results obtained from its use. GI Dynamics will not be liable for any damage, loss or expense directly or indirectly arising from the use of the device, including lost business, revenue, profits, data or goodwill. GI Dynamics assumes no liability with respect to devices that are altered, reused, reprocessed, or re-sterilized.

The end-user of the device assumes all responsibility for complying with all applicable laws, rules and regulations of all governmental authorities regarding its use.

Trademarks

EndoBarrier is a registered trademark of GI Dynamics, Inc.

For patent information, please visit:
www.gidynamics.com/patents.php

CE
0120

Not available for sale or distribution within the United States
For Export Only

MedPass International Ltd.
Windsor House
Bretforton
Evesham, Worcestershire
WR11 7JJ, United Kingdom
Tel: +44 145 2619 222

GI Dynamics, Inc.
25 Hartwell Avenue
Lexington, MA 02421 USA
Tel: +31.13.5479337

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EndoBarrier je registrovaná ochranná známka společnosti GI Dynamics, Inc.